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Hersteller

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Vertrieb

Manufacturer:
German Special Alloys GmbH | Carl-Friedrich-Benz-Straße 1b | 47877 Willich | Germany

Distributor:
Ador Dental GmbH | Zum Jägerhof 2 | 40724 Hilden | Germany

CE 0044

Instructions for Use metal-ceramic alloy

ADORBOND® BC Blank

Indications

ADORBOND® BC Blank is a dental metal-ceramic alloy based on cobalt. **ADORBOND® BC Blank** is free from nickel, cadmium, beryllium and lead and fulfils the standards of EN ISO 22674 type 4 for appliances with thin sections that are subject to very high forces, e.g. removable partial dentures, clasps, thin veneered crowns, wide-span bridges or bridges with small cross sections, bars, attachments and implant retained superstructures. This product is a milling blank intended for machining using suitable CAD/CAM milling machines (follow the manufacturer's instructions).

Contraindications

Patients with a confirmed allergy or sensitivity to a component of the alloy.
Applications not listed as indications.

Target patient groups

Patients with local tooth defects caused by dental caries or those with incomplete dentition can be treated with dental alloys for partial or full dentures. The selection of patients is the responsibility of the dentist.

Composition w_i

Co	%	61,7
Cr	%	27,8
W	%	8,5
Si	%	1,6
Nb, Mn, Fe, N	%	< 1

Technical data (standard values, milled/burned condition)

Density ρ	$\text{g} \cdot \text{cm}^{-3}$	8,4
Vickers hardness	HV 10	290
Linear CTE α (25 - 500 °C)	$10^{-6} \cdot \text{K}^{-1}$	14,1
Linear CTE α (20 - 600 °C)	$10^{-6} \cdot \text{K}^{-1}$	14,4
Highest recommended firing temperature $T_{\text{Br,max}}$	°C	980
0,2-% Yield strength $R_{p0,2}$	MPa	400
Modulus of elasticity E	GPa	195
Elongation at break A_5	%	8
Tensile strength R_m	MPa	600

Recommendations for Use

Design

The design is carried out with suitable CAD software in consideration of dental rules. Avoid wall thicknesses of less than 0,35 mm for the CAD model. Increase the wall thickness at critical points. Design the connectors as strong and high as possible (height: at least 3,5 mm, width: at least 2,5 mm).

Milling

Use suitable milling tools and cutting data according to the manufacturer's instructions for the dental milling machine.

Firing of Ceramics

Use commercially available dental ceramics for cobalt based metal alloys with a suitable linear thermal expansion coefficient. Please follow the associated work instructions and cooling schemes given by the ceramic manufacturer. After milling:

1. Separation of connectors and finishing of the object. Carbide cutters are recommended.
2. Sand blast the surface by use of a pencil-blaster with aluminium oxide 100 μm or 250 μm .
3. Ultrasonically clean the frame in distilled water or degrease with ethyl acetate.
4. The oxide firing is optional, to be done at about 960 °C under vacuum for 5 minutes. Always remove the oxide layer after oxide firing by sand blasting with aluminium oxide and degrease again.
Note: A clean surface is best to avoid bubbles in ceramics.
5. The opaque is applied on the surface by a first thin wash firing and a second evenly covering opaque layer.
Before firing always let the opaque dry for 5-10 minutes at 600 °C.
6. Firing and cooling should be carried out in accordance to the ceramic manufacturer's instructions.
7. After every firing step (dentine bake, build-up and glazing) cooling phase until ca. 750 °C.

Finishing

After firing of the ceramic, polish the frame with suitable grinding and polishing instruments for dental alloys up to high gloss.

Soldering and Welding

Soldering before firing of the frame can be carried out with commercially available solders und high temperature flux. The width of the solder gap should be 0,05-0,2 mm. For welding with laser use suitable commercially available metal welding wires.

Manufacturer:
German Special Alloys GmbH | Carl-Friedrich-Benz-Straße 1b | 47877 Willich | Germany

CE 0044

Distributor:
Ador Dental GmbH | Zum Jägerhof 2 | 40724 Hilden | Germany

Safety Note

Metal dust is harmful to health and can cause irritation if inhaled or if it comes into contact with the skin. Avoid creating dust and inhaling it! When performing finishing and sandblasting work, use a dust extraction system and wear safety goggles, protective gloves, a respirator (FFP3 – in accordance with DIN EN 149), and protective clothing. After working with the restorations, thoroughly clean the areas that have come into contact with the skin, paying particular attention to the hands. Cobalt is classified as a “CMR substance.” If ingested, cobalt is suspected of being carcinogenic, mutagenic, and toxic to reproduction. In isolated cases, electrochemically induced local irritations have been reported. Galvanic effects may occur when different alloy groups are used. Hypersensitivity to components of the alloy must be taken into account. If intolerance to individual elements of this alloy is suspected, it must not be used.

Other Notes

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the member state in which the user and/or the patient is established.

Warranty

These application recommendations are based on own experiments and experiences and can therefore only be regarded as guidelines. The dentist or dental technician is responsible for the correct processing of this alloy.

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German Special Alloys GmbH | Carl-Friedrich-Benz-Straße 1b | 47877 Willich | Germania
Ador Dental GmbH | Zum Jägerhof 2 | 40724 Hilden | Germania