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Hersteller

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Vertrieb

Manufacturer:
German Special Alloys GmbH | Carl-Friedrich-Benz-Straße 1b | 47877 Willich | Germany

Distributor:
Ador Dental GmbH | Zum Jägerhof 2 | 40724 Hilden | Germany

CE 0044

Instructions for Use metal-ceramic alloy

ADORBOND® CC Plus Powder 10 µm - 30 µm

Indications

ADORBOND® CC Plus Powder 10 µm - 30 µm is a dental metal-ceramic alloy based on cobalt. **ADORBOND® CC Plus Powder 10 µm - 30 µm** is free from nickel, cadmium, beryllium and lead and fulfils the standards of EN ISO 22674 type 5 for appliances in which parts require the combination of high stiffness and strength, e.g. thin removable partial dentures, parts with thin cross-sections and clasps. **ADORBOND® CC Plus Powder 10 µm - 30 µm** can be used as a partial denture alloy as well as a metal-ceramic alloy. This product is a metal powder intended for additive manufacturing using suitable laser sintering systems (follow the manufacturer's instructions).

Contraindications

Patients with a confirmed allergy or sensitivity to a component of the alloy.
Applications not listed as indications.

Target patient groups

Patients with local tooth defects caused by dental caries or those with incomplete dentition can be treated with dental alloys for partial or full dentures. The selection of patients is the responsibility of the dentist.

Composition w_i

Co	%	63,6
Cr	%	24,8
W	%	5,5
Mo	%	5,0
Si	%	1,1
Fe, Mn, Nb	%	< 1

Technical data (standard values, milled/burned condition)

Particle size d _i	µm	10-30
Density ρ	g · cm ⁻³	8,2
Vickers hardness	HV 10	520
Linear thermal expansion coefficient α 25 - 500 °C	10 ⁻⁶ ·K ⁻¹	14,3
Linear thermal expansion coefficient α 20 - 600 °C	10 ⁻⁶ ·K ⁻¹	14,5
Melting range TS - TL	°C	1360-1420
Highest recommended firing temperature T _{Br,max}	°C	980
0,2-% Yield strength R _{p 0,2}	MPa	1150
Modulus of elasticity E	GPa	275
Elongation at break A ₅	%	5
Tensile strength R _m	MPa	1350

Recommendations for Use

Design

The design is carried out with suitable CAD software in consideration of dental rules. Avoid wall thicknesses of less than 0,35 mm for the CAD model. Increase the wall thickness at critical points. Design the connectors as strong and high as possible (height: at least 3,5 mm, width: at least 2,5 mm).

SLM-Process

Avoid dust production while opening the powder container, transportation and filling into the SLM system! Use a suitable SLM system for the particle distribution of **ADORBOND® CC Plus Powder 10 µm - 30 µm** and carry out the SLM-Process according to the manufacturer's instructions for the SLM system. Non melted powder can be reused but should be sieved with a powder sieve (80 µm) or ultrasonic sieve (55 µm) first. Melted powder cannot be reused by e.g. casting.

Shot Peening

Blast the production platform with spherical ceramic material with Ø 0,125-0,250 mm and 2,5-3,5 bar until a uniform saturation of the surface is visibly achieved.

Stress Relief Heat Treatment

The heat treatment needs to be carried out in a special furnace under argon atmosphere. Place the production platform in an inert gas box with an argon flow of 1-2 L/min and put the box into the cold furnace. Then:

1. Within 60 min heat up to 850 °C and hold the temperature for 60 min.
2. Let the platform cool down. Open furnace door earliest at 600 °C or lower temperature, remove the inert gas box and stop the argon flow.

Manufacturer:
German Special Alloys GmbH | Carl-Friedrich-Benz-Straße 1b | 47877 Willich | Germany

Distributor:
Ador Dental GmbH | Zum Jägerhof 2 | 40724 Hilden | Germany

C € 0044

Firing of Ceramics

Use commercially available dental ceramics for cobalt based metal alloys with a suitable linear thermal expansion coefficient. Please follow the associated work instructions and cooling schemes given by the ceramic manufacturer. After the heat treatment and cooling down:

1. Remove restorations from the plate using a band saw, rotating instruments or pliers. Use pliers to remove the rest of the supports. Carbide cutters are recommended for finishing of the object.
2. Sand blast the surface by use of a pencil-blaster with aluminium oxide 100 µm or 250 µm.
3. Ultrasonically clean the frame in distilled water or degrease with ethyl acetate.
4. The oxide firing is optional, to be done at about 960 °C under vacuum for 5 minutes. Always remove the oxide layer after oxide firing by sand blasting with aluminium oxide and degrease again.
5. Note: A clean surface is best to avoid bubbles in ceramics.
6. The opaque is applied on the surface by a first thin wash firing and a second evenly covering opaque layer.
7. Before firing always let the opaque dry for 5-10 minutes at 600 °C.
8. Firing and cooling should be carried out in accordance to the ceramic manufacturer's instructions.
9. After every firing step (dentine bake, build-up and glazing) cooling phase until ca. 750 °C.

Finishing

After firing of the ceramic, polish the frame with suitable grinding and polishing instruments for dental alloys up to high gloss.

Soldering and Welding

Soldering before firing of the frame can be carried out with commercially available solders und high temperature flux. The width of the solder gap should be 0,05-0,2 mm. For welding with a laser use suitable commercially available metal welding wires.

Safety Note

Metal dust is harmful to health and can cause irritation if inhaled or if it comes into contact with the skin. Avoid creating dust and inhaling it! When handling powder, performing finishing work, or sandblasting, use a dust extraction system and wear safety goggles, protective gloves, a respirator (FFP3 – in accordance with DIN EN 149), and protective clothing. After working with powder or restorations, thoroughly clean areas that have come into contact with the skin, paying particular attention to the hands. Cobalt is classified as a “CMR substance.” If ingested, cobalt is suspected of being carcinogenic, mutagenic, and toxic to reproduction. In isolated cases, electrochemically induced local irritations have been reported. Galvanic effects may occur when different alloy groups are used. Hypersensitivity to components of the alloy must be taken into account. If intolerance to individual elements of this alloy is suspected, it must not be used.

Other Notes

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the member state in which the user and/or the patient is established.

Warranty

These application recommendations are based on own experiments and experiences and can therefore only be regarded as guidelines. The dentist or dental technician is responsible for the correct processing of this alloy.

Batch number

Follow the
instructions for use

Manufacturer

Not for reuse

Note
medical device

UDI labeling

Distributor

Manufacturer:
German Special Alloys GmbH | Carl-Friedrich-Benz-Straße 1b | 47877 Willich | Germany

Distributor:
Ador Dental GmbH | Zum Jägerhof 2 | 40724 Hilden | Germany

CE 0044

Instructions for Use partial denture alloy

ADORBOND® CC Plus Powder 10 µm – 30 µm

Indications

ADORBOND® CC Plus Powder 10 µm - 30 µm is a partial denture alloy based on cobalt. ADORBOND® CC Plus Powder 10 µm - 30 µm is free from nickel, cadmium, beryllium and lead and fulfils the standards of EN ISO 22674 type 5 for appliances in which parts require the combination of high stiffness and strength, e.g. thin removable partial dentures, parts with thin cross-sections and clasps ADORBOND® CC Plus Powder 10 µm - 30 µm can be used as a partial denture alloy as well as a metal-ceramic alloy. This product is a metal powder intended for additive manufacturing using suitable laser sintering systems (follow the manufacturer's instructions).

Contraindications

Patients with a confirmed allergy or sensitivity to a component of the alloy.

Applications not listed as indications.

Target patient groups

Patients with local tooth defects caused by dental caries or those with incomplete dentition can be treated with dental alloys for partial or full dentures. The selection of patients is the responsibility of the dentist.

Composition w_i

Co	%	63,6
Cr	%	24,8
W	%	5,5
Mo	%	5,0
Si	%	1,1
Fe, Mn, Nb	%	< 1

Technical data (standard values, milled/burned condition)

Particle size d _i	µm	10-30
Density ρ	g · cm ⁻³	8,2
Vickers hardness	HV 10	520
Melting range T _S - T _L	°C	1360-1420
0,2-% Yield strength R _{p0,2}	MPa	1150
Modulus of elasticity E	GPa	275
Elongation at break A ₅	%	5,0
Tensile strength R _m	Mpa	1350

Recommendations for Use

Digital design modelling

The design is carried out with suitable CAD software in consideration of dental rules. Wall thickness after finishing of the object: at least 0,5 mm. Base thickness of standard bases: 0,75 – 0,9 mm and 1,0 – 1,2 mm for smaller bands or skeletonised frames. Thickness of lower jaw brackets: 1,8 – 2,0 mm x 4,0 x 4,2 mm. Thickness of transition area from minor connectors to clasps: 1,5 x 2,0 mm. Thickness of clasp tip: 1,2 – 1,5 mm. Round off the areas of transition from minor connectors to clasps and the tips of clasps.

SLM-Process

Avoid dust production while opening the powder container, transportation and filling into the SLM system! Use a suitable SLM system for the particle distribution of Adorbond CC Plus Powder 10-30 µm and carry out the SLM-Process according to the manufacturer's instructions for the SLM system. Non melted powder can be reused but should be sieved with a powder sieve (80 µm) or ultrasonic sieve (55 µm) first. Melted powder cannot be reused by e.g. casting.

Shot Peening

Blast the production platform with spherical ceramic material with Ø 0,125-0,250 mm and 2,5-3,5 bar until a uniform saturation of the surface is visibly achieved.

Stress Relief Heat Treatment

The heat treatment needs to be carried out in a special furnace under argon atmosphere. Place the production platform in an inert gas box with an argon flow of 1-2 L/min and put the box into the cold furnace. Then:

1. Within 90 min heat up to 1100 °C and hold the temperature for 60 min.
2. Let the platform cool down. Open furnace door earliest at 600 °C or lower temperature, remove the inert gas box and stop the argon flow.

Manufacturer:
German Special Alloys GmbH | Carl-Friedrich-Benz-Straße 1b | 47877 Willich | Germany

Distributor:
Ador Dental GmbH | Zum Jägerhof 2 | 40724 Hilden | Germany

CE 0044

Treatment of the restoration

After the heat treatment and cooling down:

1. Remove restorations from the plate using a band saw, rotating instruments or pliers. Use pliers to remove the rest of the supports. Carbide cutters are recommended for finishing of the object.
2. Sand blast the surface by use of a pencil-blaster with aluminium oxide 100 µm or 250 µm.
3. Electrolytic polishing with commercial electrolytes in dental polishing units. Cover clamps and fitting parts with covering varnish while polishing.

Finishing

After finishing and fitting smooth the frame with a rubber polisher.

Soldering and Welding

Soldering before firing of the frame can be carried out with commercially available solders und high temperature flux. The width of the solder gap should be 0,05-0,2 mm. For welding with a laser use suitable commercially available metal welding wires.

Safety Note

Metal dust is harmful to health and can cause irritation if inhaled or if it comes into contact with the skin. Avoid creating dust and inhaling it! When handling powder, performing finishing work, or sandblasting, use a dust extraction system and wear safety goggles, protective gloves, a respirator (FFP3 – in accordance with DIN EN 149), and protective clothing. After working with powder or restorations, thoroughly clean areas that have come into contact with the skin, paying particular attention to the hands. Cobalt is classified as a “CMR substance.” If ingested, cobalt is suspected of being carcinogenic, mutagenic, and toxic to reproduction. In isolated cases, electrochemically induced local irritations have been reported. Galvanic effects may occur when different alloy groups are used. Hypersensitivity to components of the alloy must be taken into account. If intolerance to individual elements of this alloy is suspected, it must not be used.

Other Notes

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the member state in which the user and/or the patient is established.

Warranty

These application recommendations are based on own experiments and experiences and can therefore only be regarded as guidelines. The dentist or dental technician is responsible for the correct processing of this alloy.

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German Special Alloys GmbH | Carl-Friedrich-Benz-Straße 1b | 47877 Willich | Germania
Ador Dental GmbH | Zum Jägerhof 2 | 40724 Hilden | Germania